



**PreserveTech™**  
Stays Cleaner • Lasts Longer

# Secure Lock Raised Toilet Seat

The new PreserveTech™ Secure Lock Raised Toilet Seat features three locking points on the front and sides for added security and stability. This toilet seat adds an extra 5" of height to the toilet and includes padded armrests that can be adjusted or removed.

Most non-permanent raised toilet seats only offer one or two locking points of security resulting in instability when the end user places more weight on one side of the unit. The seat fits on most elongated and standard toilets and installs without tools.

Drive PreserveTech™ products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

\*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.



**HOW TO CHOOSE A RAISED TOILET SEAT**



*Scan me*

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
RTL12C003-WH	PreserveTech™ Secure Lock Raised Toilet Seat	822383005669/ea / 50822383005664/cv	1/ea, 1/cv

**PRODUCT SPECS**

Product Dimensions	20.32" (W) x 17.44" (L) x 11.06" (H)
Product Weight	5.29 lb.
Color	White
Weight Capacity	350 lb.
Warranty	Limited Lifetime**

**BOX SPECS**

Retail Box Dimensions	20.87" (W) x 7.48" (L) x 18.11" (H)
Retail Box Weight	5.95 lb.

\*\*Refer to instruction guide or website for warranty details.