

Aquasense Adjustable Transfer Bench

HCPCS: E0247,E0247

- Both lightweight and durable the rust-resistant aluminum frame supports up to 250 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- A removable push-button armrest supports transferring in and out of the bathtub.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 18" to 22".
- Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
- Push button, snap together design requires no tools and assembles in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.

PRODUCTS

Description

770-406	AQUASENSE TRANSFER BENCH, 1EA, 2CV, Ret
---------	-----------------------------------------

SPECIFICATIONS

Description

Clearance	14.5 to 18.5 (36.8 cm to 47 cm)
Height Adjustment	4 in 0.47 (1.2 cm) increments
Overall Height	30 to 34 (76.5 to 86.7 cm)
Overall Width	26 (66 cm)
Retail Box Dimensions	20.9 x 17.7 x 19.9 53 cm x 45 cm x 50.5 cm
Seat To Floor Height	18 to 22 (45.7 to 55.9 cm)
Waranty	Limited Lifetime
Weight Capacity	250 lb (113 kg)