

## PreserveTech Universal Raised Toilet Seat



HCPCS: E0244



FIGURE A



FIGURE B



FIGURE C

- Accommodates three different seat heights (2", 4" and 6") by maneuvering just two components.
- Pivoting armrests can be raised out of the way or removed depending on your needs or bathroom environment.
- Fits comfortably on most standard and elongated toilets.
- Includes a lid with recessed grip for better hygiene and aesthetic and effortless opening.
- Contoured to provide a safe and comfortable elevation.
- Made from PP that is both lightweight and durable to support up to 300 lbs.
- Drive PreserveTech products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.
- \*Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

### PRODUCTS

PRODUCTS	Description
RTL12C002-WH	PreserveTech Universal RTS, 1EA, 1/CV

### Description

### SPECIFICATIONS

SPECIFICATIONS	Description
Seat Depth	16"
Seat Width	19"
Warranty	Limited Lifetime
Weight Capacity	300 lbs
Width Between Arms	18.5"
Material Composition	Polypropylene