

Folding Shower Chair

HCPCS: E0245



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Product available with back (main image) and without (Figure A)
- Conveniently folds flat (Figure B)
- Blow-molded bench provides comfort and strength
- Drainage holes in seat reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Angled legs provide additional stability
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
12486	Bath Bench,Folding,1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	6.6 lbs. 7 lbs.
Outside Legs	20.5" (W) x 19.25" (D) 17.75" (W) x 14" (D)
Seat Dimensions	19.75" (W) x 12" (D) 18" (W) x 11" (D)
Seat Height	17.5" - 18.5" 16.25" - 18.25"
Warranty	Limited Lifetime
Weight	4.5 lbs. 5 lbs.
Weight Capacity	300 lbs

ACCESSORIES	Description
12486T	Replacement Tips for 12486 4/Bx, BX