

Combination Transfer Bench/Commode

HCPCS: E0247



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Tool-free assembly of back, legs and arm (Figure A - C)
- Combines transfer bench and commode into one
- Durable blow-molded plastic bench and backrest
- Aluminum A-frame construction is lightweight and provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- No exposed hardware to injure patient
- Extra-large suction cups provide added safety
- New pinch-free cover allows for push pins to be depressed without pinching fingers (Figure D)

PRODUCTS

Description

| | |
|------------|--------------------------------|
| 12011KDC-1 | Tranf Bnch w/Commode, 3pc,1/cs |
|------------|--------------------------------|

SPECIFICATIONS

Description

| | |
|------------------------|-------------------------|
| Carton Shipping Weight | 13 lbs. 24 lbs. |
| Outside Legs | 31.5" (W) x 22.25" (D) |
| Seat Dimensions | 26.25" (W) x 16.25" (D) |
| Seat Height | 18" - 22" |
| Warranty | Limited Lifetime |
| Weight | 11 lbs. |
| Weight Capacity | 400 lbs |

ACCESSORIES

Description

| | |
|---------|---------------------------------|
| 12011SL | SuctionCups,Lg.12011,1pr/bx, BX |
| 12011SS | SuctionCups,Sm.12011,1pr/bx, BX |