

Folding Universal Sliding Transfer Bench

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Transfer Bench slides left and right along frame for easy entry into and out of tub
- Seat height adjusts in $\frac{1}{2}$ " increments
- Folds flat for transport and storage (Figure A)
- Comes standard with removable soap dish (Figure B)
- Comes with safety net to catch dropped personal items and prevent bending
- Comes with seat belt for added safety
- Retail packaged (Figure D)
- Made using recycled materials

PRODUCTS

	Description
RTL12075	Sliding Shower Chair w/Soap Box,1EA,2/CS

SPECIFICATIONS

	Description
Carton Shipping Weight	38 lbs.
Height	21" - 25"
Seat Dimensions	22.5" (W) x 18.5" (D)
Warranty	Limited Lifetime
Weight	16 lbs.
Weight Capacity	300 lbs

ACCESSORIES

	Description
12075-A	Arm,RTL12075,1ea, EA