

## Full-Electric Bariatric Bed, 48"

HCP/CS: E0302/E0304/E1399



FIGURE A

- Quiet, smooth operation
- Heavy-duty frame ensures strength and patient safety
- All steel construction
- Larger sleep surface than conventional bed
- Split-pan design with removable bed ends is easy to set up
- No tools necessary for set up
- Hand control (included) provides for multiple bed positioning for patients
- Emergency manual crank included
- Low air loss/alternating pressure mattress available
- "T" rails available (Sold separately) (Figure A)
- Weight Capacity 750 lbs

### PRODUCTS

	Description
15302	Bariatric Bed,FullElectric,48"1/cs
15302BV-1HR	BariatricBed,FullElectric,48",1/cs
15302BV-PKG	BariatricBed,FullElectric,48",1/cs

### SPECIFICATIONS

	Description
Bed Height	18.25" - 26.25"
Carton Shipping Weight	343 lbs.
Unit Dimensions	48" (W) x 88" (L)   Overall Sleep Surface 48" W x 80" L
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components
Weight	312 lbs.
Weight Capacity	Total Weight Capacity 750 lbs., Maximum Patient Weight 650 lbs.
Bed Min height	18.25"
Bed Max height	26.25"

### PARTS

	Description
15350H	Head Board for 15302,1/cs, CS
15309N	HeadFrame,15302,BariatricBed,1/cs, CS
15320N	FootFrame,15302,BariatricBed,1/cs, CS
15350F	Foot Board for 15302,1/cs, CS
15302P1005	CompositeHeadPanels,15302,1/pr/cs, CS
15302P1006	CompositeFootPanel,15302,1pr/ea, EA
15302P1051	PanelFrame,HeadBoard,15302,1/ea, EA
15302P1052	PanelFrame,FootBoard,15302,1/ea, EA

PARTS	Description
15032MHN	HeadMotor15302BariatricBed,NewStyle,1/ea, EA
15302JB-4PORT	4 Port Junction Box for 15302/15303,1ea, EA
15300P1016	Cotter Pin for 15300/2/3 Beds, 1/ea, EA
15300CCN	Motor Connection Cord for 15300, 1/ea, BX
15302HL	HiLoMotor15302/15303,1ea, EA
15300FSR	FrameSupportRails,15302,15303,1/ea, EA
15300CL	Caster w/Lock,15300,15302,15303,1/ea, EA
15300P1015L	Adjust.Leg,Baria.Bed Headboard,1/ea, EA
15300PH1014	Insert,Head Tube,Bariat.Bed,1/ea, EA
15300EC	Caster Socket7/8 Dia shaft 1/ea, EA
15300MF	Motor,Foot,Baritaric,1/ea, EA
15005P-1028	Quick Release for Bed Motor,1/ea, EA