

Full-Electric Bariatric Bed, 48"

HCPCS: E0302




FIGURE A

- Quiet, smooth operation
- Heavy-duty frame ensures strength and patient safety
- All steel construction
- Larger sleep surface than conventional bed
- Split-pan design with removable bed ends is easy to set up
- No tools necessary for set up
- Hand control (included) provides for multiple bed positioning for patients
- Emergency manual crank included
- Low air loss/alternating pressure mattress available
- "T" rails available (Sold separately) (Figure A)
- Weight Capacity 750 lbs

PRODUCTS

PRODUCTS	Description
15302	Bariatric Bed,FullElectric,48"1/cs
15302BV-1HR	BariatricBed,FullElectric,48",1/cs
15302BV-PKG	BariatricBed,FullElectric,48",1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Bed Height	18.25" - 26.25"
Carton Shipping Weight	343 lbs.
Unit Dimensions	48" (W) x 88" (L) Overall Sleep Surface 48" W x 80" L
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components
Weight	312 lbs.
Weight Capacity	Total Weight Capacity 750 lbs., Maximum Patient Weight 650 lbs.
Bed Min height	18.25"
Bed Max height	26.25"

PARTS

PARTS	Description
15350H	Head Board for 15302,1/cs, CS
15309N	HeadFrame,15302,BariatricBed,1/cs, CS
15320N	FootFrame,15302,BariatricBed,1/cs, CS
15350F	Foot Board for 15302,1/cs, CS
15302P1005	CompositeHeadPanels,15302,1/pr/cs, CS
15302P1006	CompositeFootPanel,15302,1pr/ea, EA
15302P1051	PanelFrame,HeadBoard,15302,1/ea, EA
15302P1052	PanelFrame,FootBoard,15302,1/ea, EA

PARTS	Description
15032MHN	HeadMotor15302BariatricBed,NewStyle,1/ea, EA
15302JB-4PORT	4 Port Junction Box for 15302/15303,1ea, EA
15300P1016	Cotter Pin for 15300/2/3 Beds, 1/ea, EA
15300CCN	Motor Connection Cord for 15300, 1/ea, BX
15302HL	HiLoMotor15302/15303,1ea, EA
15300FSR	FrameSupportRails,15302,15303,1/ea, EA
15300CL	Caster w/Lock,15300,15302,15303,1/ea, EA
15300P1015L	Adjust.Leg,Baria.Bed Headboard,1/ea, EA
15300PH1014	Insert,Head Tube,Bariat.Bed,1/ea, EA
15300EC	Caster Socket7/8 Dia shaft 1/ea, EA
15300MF	Motor,Foot,Baritaric,1/ea, EA
15005P-1028	Quick Release for Bed Motor,1/ea, EA