

## Tool-Free Adjustable Length Home-Style Bed Rail

HCPCS: E0310





FIGURE A




FIGURE B



FIGURE C

- Rail length can be adjusted without tools to 37", 43" or 57" long with simple push buttons (Figure A , B)
- Crossbars are easy-to-install and remove Adjusts between 36" and 72" Accommodates a KING SIZE BED Comes standard
- Fits all home-style beds (Figure C)
- Easily adjusts up or down with spring-loaded release
- Attractive brown-vein finish
- Constructed of 1" steel

### PRODUCTS

#### Description

16500BV	Bed Rail, Home Style, Brown Vein, 1pr/cs
---------	--

### SPECIFICATIONS

#### Description

Crossbar Adjustment	36"-72"
Rail Length Adjustment	37", 43", or 57"
Warranty	Limited Lifetime

### ACCESSORIES

#### Description

16500CBV	Cross Barsw/extensions,16500BV, 1pr/cs, CS
----------	--