

Multi-Height Manual Bed

HCPCS: E0293




FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Convenient fold down hand cranks at the foot end for Foot, Head and overall bed height adjustment
- Easy to set up, bed ships in two cartons, shaft stores under frame when unassembled
- Head/foot adjustment provide anatomically correct sleep surface
- Reinforced channel frame construction offers superior strength, resists twisting and reduced weight
- Long lasting Zinc-coated spring deck for Pt comfort
- Heavy Duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- Satisfies applicable FDA entrapment guidelines
- Patient weight Capacity 350 lbs. (158kg). Total Weight Capacity 450 lbs. (204kg)
- Half rails and full rails available in packages (PKG) or sold separately

PRODUCTS	Description
15003P	Manual Bed,1/cv
15003P-FR	Manual Bed w/Full Rails Pkg
15003P-FRPKG	Manual Bed FRail IMatt Pkg
15003P-FRPKG-FM	Manual Bed FRail FMatt Pkg
15003P-FRPKG-TS	Manual Bed FRail TSMatt Pkg
15003P-HR	Manual Bed w/Half Rails Pkg
15003P-HRPKG	Manual Bed HRail IMatt Pkg
15003P-HRPKG-TS	Manual Bed HRail TSMatt Pkg

SPECIFICATIONS	Description
Bed Height	14.96" (380mm) - 22.83" (580mm) deck height
Suitable Room Size	36" (.914m) (W) x 88" (2.24m) (L)
Unit Dimensions	36" (.914m) (W) x 88" (2.24m) (L)
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components
Weight	156 lbs. (71kg)
Weight Capacity	Pt. Weight 350 lbs. (158kg) Total Weight 450 lbs (204kg) includes weight of occupant, all accessories such as mattress, side/assist rails, trapeze, any personal affects, etc.