

## Semi-Electric Bed (Single Crank)

HCPCS: E0295




FIGURE A



FIGURE B

- Quiet, smooth operation, Head/Foot section adjusts electronically, Height adjusts manually
- Easy-to-transport and set up, Bed ships in two cartons
- Shaft stores under frame when unassembled
- Channel frame construction provides superior strength and reduced weight
- Reinforced frame resists twisting and bending.
- Long Lasting Zinc-coated spring deck for Pt comfort
- Heavy Duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- Clip-on pendant mounts on bed rail, operate Head and Foot section simultaneously
- Half rails and full rails available in packages (PKG) or sold separately
- Satisfies applicable FDA entrapment guidelines
- Extension kit available (model #15005EXTKIT-L) extends length from 80" to 84"
- Patient weight Capacity 350 lbs. (158kg). Total Weight Capacity 450 lbs. (204kg)
- Motor is UL Approved, Bed is manufactured in ISO9001 & ISO13485 facility

### PRODUCTS

#### Description

15004P	SemiElec Bed,1/cv
15004P-FR	SemiElec Bed w/Full Rails Pkg
15004P-FRPKG	SemiElec Bed,FRail IMatt Pkg
15004P-FRPKG-FM	SemiElec Bed,FRail FMatt Pkg
15004P-HR	SemiElec Bed w/Half Rails Pkg
15004P-HRPKG	SemiElec Bed HRail IMatt Pkg
15004P-HRPKG-TS	SemiElec Bed HRail TSMatt Pkg

### SPECIFICATIONS

#### Description

Bed Height	14.96" (380mm) - 22.83" (580mm) deck height
Carton Shipping Weight	176 lbs. (80kg)
Unit Dimensions	36" (.914m) (W) x 88" (2.24m) (L)   Overall Sleep Surface 36" W x 80" L
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, motor, pendant, 1 Year on All Other Parts and Component
Weight	161 lbs. (73kg)
Weight Capacity	Pt. Weight 350 lbs. (158kg) Total Weight 450 lbs (204kg) includes weight of occupant, all accessories such as mattress, side/assist rails, trapeze, any personal affects, etc.