

Deluxe Bariatric Drop-Arm Commode

HCPCS: E0168



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy-to-maintain
- Large, durable snap-on seat
- Padded armrests provide comfort and stability
- Pail easily slides in from the front of the commode
- Comes complete with 12 qt commode bucket and cover that removes from front
- Weight Capacity 1000 lbs
- Made using recycled materials

PRODUCTS	Description
11135-1	Commode,Oversized DropArm,1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	26 lbs.
Height	17.5" - 22"
Seat Dimensions	23.25" (W) x 18.5" (D)
Warranty	Limited Lifetime
Weight	21 lbs.
Weight Capacity	1000 lbs
Width	25.75" between arms 26.25" Outside

PARTS	Description
RTL10389G	Utility ReplaceTips,1",Grey,RTL,1PR,8/CS, CS
11135B	Back Support Only,11135,1/ea, EA
11120B	Commode Bucket for 11120, 12qt, 1/ea, EA
11135L	Leg for 11135 1Z Serial#,1/ea, EA
11135A	Arm Only for 11135,1/ea, EA
111351Z-S	Seat Only 1Z, 11135,1/ea, EA
11135BRKN	Pail Bracket Set,New,11130/11135,1pair, PR
RTL10389GB	Utility Replac Tip, 1",Gry, 1PC/BL,10/, EA
12002PB	Push Button for Bath Bench,Walker,1/ea, EA