

Deluxe Steel Drop-Arm Commode with Padded Seat

HCPCS: E0165



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Easy-to-assemble without tools in 4 easy steps (Figure A - D)
- Padded open front vinyl toilet seat
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode
- Easy-to-clean gray powder-coated steel finish
- Padded armrests and back provide extra comfort and support
- Comes with 12 qt commode bucket with carry handle and splash shield
- Made using recycled materials

PRODUCTS

Description

| | |
|-------------|-------------------------------------|
| 11125PSKD-1 | Commode,DropArm,Padded Seat,KD,1/cs |
|-------------|-------------------------------------|

SPECIFICATIONS

Description

| | |
|------------------------|-----------------------------------|
| Carton Shipping Weight | 23 lbs. |
| Height | 19" - 23" |
| Seat Dimensions | 15.25" (W) x 16.75" (D) |
| Warranty | Limited Lifetime |
| Weight | 21 lbs. |
| Weight Capacity | 300 lbs |
| Width | 19" between arms 23.75" Outside |

PARTS

Description

| | |
|--------------|--|
| 11125PSKD-AL | LeftArmWithReleaseLever11125PSKD,1ea, EA |
| 11125L | Replacement leg, 11125,1/ea, EA |
| 12002PB | Push Button for Bath Bench,Walker,1/ea, EA |
| 11125PSKD-B | Back for 11125PSKD, 1/ea, EA |