

Respiratory Rollator Accessory Pack



HCPCS:



FIGURE A



FIGURE B

- Removable oxygen carry bag
- Multipurpose mount for cell phones, smart phones, and pulse oximeters
- Tubing and wire system to keep oxygen tubing and pulse oximeter wire probes from getting snagged
- Can fit up to two oxygen cylinders (ML6, M4, M9, M6) or one portable oxygen concentrator

PRODUCTS

RTL10555RD-RESP

Description

I-walker w/Resp Kit, Red

SPECIFICATIONS

Description

Dimensions 16" (W) x 6" (L) x 8.5" (H)

Handle Height 33.75" - 38.75"

Seat Dimensions 18" (W) x 9.75" (D)

Seat Height 20"

Unit Dimensions 27" (L) x 25.625" (W)

Warranty 1 Year Limited | Limited Lifetime | 5 Years

Weight 1.4 lbs. | 17 lbs.

Weight Capacity 300 lbs

Weight Capacity (13046 | 13046-HD) 300 lbs.

specifications I-Walker rollator | accessory pack carry bag