

Travelite Transport Chair



HCPCS:



FIGURE A



FIGURE B




FIGURE C



FIGURE D

- The Travelite Chair in a Bag is a lightweight yet sturdy aluminum transport chair in gunmetal gray
- It is easily transportable and comes with a compact carry bag with wheels making it ideal for traveling (Figure A - B)
- Folds for storage and transportation (Figure C)
- It is ideal for use indoors and outdoors
- 18" square nylon seat comes with lift-up straps
- Comes with wheel locks 8" solid tires
- Comes with padded, flip-up desk arms and height adjustable footrests

PRODUCTS	Description
TC005GY	Travelite Chair in a Bag,1/ea

SPECIFICATIONS	Description
Back Height	17.5"
Carton Shipping Weight	23.2 lbs
Overall Dimensions	22.5" (W) (open) 11.5" (W) (closed) x 37.5" (H)
Overall Length (with riggings)	33"
Seat Dimensions	18" (W) x 18" (D)
Seat To Floor Height	19"
Warranty	Limited Lifetime
Weight	20.5 lbs
Weight Capacity	250 lbs

PARTS	Description
TC00503	Upholstery,TC005,1ea, EA
TC00504L	ArmRestAssemblyLeftCompleteTC005,1ea, EA
TC00509	Foot Rest Connector,TC005,1ea, EA
TC005012	Foot Rest Assembly,TC005,1pair, PR
TC005027	Lifting Strap for TC005GY, 1/ea, EA
TC00507	Arm Rest Release Knob,TC005,1ea, EA
TC005011	CrossFrameBlocksSetof4TC005,1package, PKG
STDS1060	Wheelchair Fork Bearing,1/ea, EA
TC005020	FrontForkTC005,1ea, EA
STDS4S3602	Rear Axle Bearing,TR4S,1pair, PR
TC005021F	FrontWheelWithBearingsTC005,1ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS3J2422	Frame Cap,Cruiser3,1/ea, EA

PARTS

Description

TC005024

Carry Bag,TC005,1/ea, EA