

## 3 Drawer Chest

HCPCS:

- High Pressure wood grained Laminate Cabinet finish, Buried wood edges with High Pressure Laminate Top, Drawer front and door feature raised real wood detail, Solid wood feet
- Available in six most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- Ventilated back allows air flow throughout cabinet interior (Figure B)
- Available in Arched and Flat Toe Kicks (Figure C)

<b>PRODUCTS</b>	<b>Description</b>
RS-CHST-03	Riverside Chest 3drwr
RS-CHST-03W	Riverside Chest 3drwr-wide