

## Hamper Stand

HCPCS:



- Supports 36 42 gallon bags
- Glides easily with 3" locking casters
- Chrome-plated steel is attractive and easy-to-maintain
- Non-slip pedal raises the white, powder-coated steel lid
- Ships knocked down

### PRODUCTS

	Description
13070	Hamper Stand, 1/cs

### SPECIFICATIONS

	Description
Base	16.25" (W) x 19.25" (D)
Carton Shipping Weight	28 lbs.
Height	38"
Lid Dimensions	19" (W) x 19.5" (D)
Warranty	Limited Lifetime
Weight	24 lbs.