

## Hydraulic, Deluxe Silver Vein Patient Lift

HCPCS: E0630



- Silver Vein steel construction
- 6 point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base

**PRODUCTS****Description**

13023SV	Patient Lift, New Style, Silver Vein, 1/cs
---------	--

**SPECIFICATIONS****Description**

Base Height	6.5" (13023SV)
Base Width	42" (open) / 22.5" (closed)
Boom Height	17.5" - 72" (13023SV)
Carton Shipping Weight	84 lbs.
Casters	5" (13023SV)
Warranty	Limited Lifetime (frame), 1 Year Limited (pump)
Weight	72 lbs.
Weight Capacity	450 lbs

**PARTS****Description**

13023SV01	Hanger Attachment Bolt 13023SV, 1 ea, EA
13023SV03	Mast Complete 13023SV, 1 ea, EA
13023SV04	Upper Mast 13023SV, 1 ea, EA
13023SV06	Pump Attachment Bolt 13023SV, 1 ea, EA
13023SV07	Mast Hinge Bolt 13023SV, 1 ea, EA
13023SV08	Push Handle Grip 13023SV, 1 PR=1 EA, PR
13023SV09	Pump 13023SV, 1 ea, EA
13023SV09-FPHH	Pump Handle Hardware, 13023SV, 1/ea, EA
13023SV10	Spreader Bar Handle Grip 13023SV, 1 ea, EA
13023SV11	Spreader Bar 13023SV, 1 ea, EA
13023SV12	Thumb Screw 13023SV, 1 ea, EA
13023SV13	Spreader Bar Receiver 13023SV, 1 ea, EA
13023SV16	Rear Caster With Lock 13023SV, 1 ea, EA
13023SV18	Front Caster W/o Lock 13023SV, 1 ea, EA
13023SV19	Base Tip Cover 13023SV, 1 ea, EA
13023SV20	Tube Plug 13023SV, 1 ea, EA

**PARTS****Description**

13023SV21	CradlePivot13023SV,1ea, EA
13023SV22	PumpHandGrip13023SV,1ea, EA
13023SV23	MastTipCover13023SV,1ea, EA
13052RL	ReleaseLever,13052,1/ea, EA