

Fish Pediatric Compressor Nebulizer

HCPCS: E0570-1,E0570



FIGURE A



FIGURE B



FIGURE C

- Child-friendly design to encourage therapy compliance
- Easy-to-replace filter
- Powerful piston pump is efficient, dependable and quiet
- Complete with disposable neb kit, mouth-piece, pediatric aerosol mask, carry bag, replacement filters and instructions
- Retail packaged

PRODUCTS

Description

18090-FS	Ped FishNeb w/Bag,DisposNebKit,1/ea,6/c
18091-FS	PedFishNebw/BagReuseDisposNeb,1/ea,6/cs

SPECIFICATIONS

Description

Electrical Requirements	120 VAC 60 Hz
Liter Flow	8 lpm
Maximum Pressure	30 psi
Method of Nebulization	Piston Compressor
Operating Pressure	10 psi
Particle Size	0.5 µm to 10 µm
Unit Dimensions	5.25" (W) x 6.75" (H) x 9.25" (D)
Warranty	Limited 5 Year
Weight	3.1 lbs

ACCESSORIES

Description

NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
18254	Reusable Neb Kit 10/case, CS
NEB KIT 700	NebKit,IncPediMask,Tub ,50/cs, CS