

## Panda Pediatric Compressor Nebulizer

HCPCS: E0570, A7003, A7005, A7015, A7013



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Child Friendly Design
- 6003 and 6005 models include disposable neb kit, pediatric mask, tubing, and a carry bag (6004 models exclude carry bag)
- Durable, piston driven system

### PRODUCTS

#### Description

MQ6003	Panda Compressor Nebulizer 1/ea, 6/cv
MQ6003R	PandaNebw/Reuse&DisposNeb,1/ea, 6/cv
MQ6004	Panda Compressor Nebw/out Bag 1/ea, 8/cv
MQ6005	Pink Panda Compressor Neb, 1/ea, 6/cv
MQ6005R	Pink Panda Neb, w/ReuseNebCup 1/e, 6/cv

### SPECIFICATIONS

#### Description

Compressor Pressure Range	30 to 36 psi
Liter Flow Range	4 to 8L/min.
Nebulizer Medication Capacity	5mL
Noise Level	55 dBA
Operating Humidity Range	10% to 95% RH
Operating Pressure Range	8 to 16 psi
Operating Temperature Range	10° to 40°C (50° to 104°F)
Particle Size	0.5 µm to 5.0 µm
Power Consumption	80 VA
Power Requirements	120V, 50-60 Hz
Storage Humidity Range	10% to 95% RH
Storage Temperature Range	-20° to 70°C (-4° to 158° F)
Warranty	Limited 5 Year
Weight	4.3 lbs

### ACCESSORIES

#### Description

NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
18254	Reusable Neb Kit 10/case, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS