

Hybrid Universal Interface

HCPCS:



- Universal design simplifies choice and fitting protocols and minimizes time spent with a patient to address issues such as mouth breathing or claustrophobia
- Adjustable nasal pillow design provides a customized fit for each patient
- No nasal bridge or forehead contact eliminates nasal bridge sores and air leaks near eyes
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Built-in chin flap gently supports the chin and keeps the mouth closed, eliminating the need for a chin strap
- Premium headgear provides five points of adjustability to provide better fit, seal and increased support

SPECIFICATIONS

Description

Warranty	90 days
----------	---------

ACCESSORIES

Description

HYB513	Hybrid Oral Cus - Med, EA
HYB515	Hybrid Oral Cus - Lg, EA
HYB521	Hybrid Nasal Pillows - Sm, EA
HYB530	Hybrid HdGr, EA
HYB523	Hybrid Nasal Pillows - Med, EA
HYB513	Hybrid Oral Cus - Med, EA
HYB515	Hybrid Oral Cus - Lg, EA
HYB521	Hybrid Nasal Pillows - Sm, EA
HYB530	Hybrid HdGr, EA
HYB523	Hybrid Nasal Pillows - Med, EA